

**TECHNICAL CARD: SPYLOG**

<b>DESCRIPTION</b>	Software for monitoring continuously analog values, according to the 21 CFR Part 11 FDA regulation
<b>FEATURES</b>	<ol style="list-style-type: none"><li>1. <i>Complete</i>: SPYLOG has been developed in order to manage the values acquired from the field and the indication of alarm events completely (processing, saving, displaying and printing)</li><li>2. <i>Capable</i>: the system manages up to 80 channels, as well as, their relative alarm and prealarm thresholds overflow</li><li>3. <i>Multilingual</i>: the product is available in Italian, English and German</li><li>4. <i>Certified</i>: the software is certified according to the "21 CFR Part 11", the FDA regulation for electronic records and their authentication through the use of the electronic signature</li></ol>
<b>FUNCTIONALITIES</b>	<ol style="list-style-type: none"><li>1. Display of the values acquired from the field (Trend and Bar-graph)</li><li>2. Indication of alarms</li><li>3. Printing of reports and graphic pages</li><li>4. Display of communication errors</li><li>5. Backup, restore and cancellation of historical data</li></ol>
<b>RESTRICTIONS</b>	The SPYLOG current version does not manage yet data exportation to other systems
<b>SYSTEM REQUIREMENTS</b>	<ol style="list-style-type: none"><li>1. Windows NT4 or Windows 2000 Professional Operating Systems</li><li>2. 32 Mbyte RAM</li><li>3. Hard Disk of 2 Gbyte</li><li>4. CD-RW Drive</li><li>5. Ethernet Card</li></ol>
<b>PRODUCT CONTENT</b>	<p>The package available in CD contains:</p> <ol style="list-style-type: none"><li>1. <b>PRODUCT</b>: SPYLOG</li><li>2. <b>USER GUIDE</b>: user and installer guides</li></ol> <p>If willing to buy also the data acquisition system along with the material previously described, the package will also include:</p> <ol style="list-style-type: none"><li>1. <b>PRODUCT</b>: iBOX MASTER and REMOT + Setup software</li><li>2. <b>USER GUIDE</b>: installation guide</li></ol>