

TECHNICAL CARD: SPYLOG

DESCRIPTION	Software for monitoring continuously analog values, according to the 21 CFR Part 11 FDA regulation
FEATURES	<ol style="list-style-type: none">1. <i>Complete</i>: SPYLOG has been developed in order to manage the values acquired from the field and the indication of alarm events completely (processing, saving, displaying and printing)2. <i>Capable</i>: the system manages up to 80 channels, as well as, their relative alarm and prealarm thresholds overflow3. <i>Multilingual</i>: the product is available in Italian, English and German4. <i>Certified</i>: the software is certified according to the "21 CFR Part 11", the FDA regulation for electronic records and their authentication through the use of the electronic signature
FUNCTIONALITIES	<ol style="list-style-type: none">1. Display of the values acquired from the field (Trend and Bar-graph)2. Indication of alarms3. Printing of reports and graphic pages4. Display of communication errors5. Backup, restore and cancellation of historical data
RESTRICTIONS	The SPYLOG current version does not manage yet data exportation to other systems
SYSTEM REQUIREMENTS	<ol style="list-style-type: none">1. Windows NT4 or Windows 2000 Professional Operating Systems2. 32 Mbyte RAM3. Hard Disk of 2 Gbyte4. CD-RW Drive5. Ethernet Card
PRODUCT CONTENT	<p>The package available in CD contains:</p> <ol style="list-style-type: none">1. PRODUCT: SPYLOG2. USER GUIDE: user and installer guides <p>If willing to buy also the data acquisition system along with the material previously described, the package will also include:</p> <ol style="list-style-type: none">1. PRODUCT: iBOX MASTER and REMOT + Setup software2. USER GUIDE: installation guide